



# PRODUCT CATALOG

**Animal  
System**



# Animal Whole-body plethysmography System

The whole-body plethysmography (WBP) system allows for the measurement of lung function and airway response in awake, freely moving animals. It avoids the effects of invasive tracheotomy and anesthesia, making the procedure easy and fast, and is suitable for long-term follow-up studies.



The animal is placed in whole-body plethysmograph, which is connected to a highly sensitive sensor. When the animal breathes, its thorax undulates and changes the volume of the body tracing chamber. The volume change is converted into electrical signal by pressure transducer and amplifier, and the respiratory curve is displayed on the computer screen after processing. Some respiratory parameters can be calculated, such as the tidal volume (TV), peak expiratory flow rate (PEF), respiratory rate, etc.

## Simple to operate

- No anesthesia, No surgery, Free moving.
- Anesthesia and surgery can affect the animal's spontaneous breathing.

# ADVANTAGES

## High throughput

- A maximum of 64 animals can be monitored simultaneously.

## Long-term tracking study

- Animals are in awake and free-moving state during measurement, this allows for long-term follow-up studies, suitable for initial drug screening


## Acquire the real data

- Non-invasive method can obtain the respiratory rate, tidal volume and minute ventilation under the real spontaneous breathing. The invasive mode is the animals applied by mechanical ventilation after anesthesia surgery, so the above parameters cannot be obtained accurately.




# Animal Whole-body plethysmography System

## Measuring parameters




- Breath Frequency
- Tidal Volume
- Respiratory Flow
- Minute Ventilation
- Inspiratory Time
- Expiratory Time
- Peak Inspiratory Flow Rate
- Peak Expiratory Flow Rate
- End Inspiratory Pause
- End Expiratory Pause
- Relaxation Time
- Penh
- Accumulated Volume
- VO2 (optional)
- VCO2(optional)
- RER(optional)
- EE(optional)
- Other Parameters

## Cough detection



# Animal Whole-body plethysmography System


## Application



## Drug aerosol administration



The system can be extended to connect with a nebulizer to administer drug to the animals, this is generally used in cough and asthma studies. WBP system can detect and record the respiratory parameters and airway hyperresponsiveness during excitation. The WBP system can detect and record the changes of respiratory indexes and airway hyper-responsiveness during the excitation process, which can be used to evaluate the broncho-constriction, etc.



## Software

### Easy Setup

Automatic calibration  
Complete calibration operations  
in 1 minute.

### Powerful

Powerful analytical capabilities, the software can perform both single-wave analysis and trend analysis, it greatly facilitates data processing

### Data&Report

Each waveform data and analyzed data can be stored and exported in the format of txt, csv, pdf, jpg

### GLP Features

- Control over User Permissions
- Audit Trails
- Electronic Signatures
- Data Integrity
- Electronic Signatures


## Optional Plethysmographs



With water bottles and food trough, available for long-term continuous monitoring



Electrophysiology plethysmograph, for combined with optogenetic technology, EEG, EMG, etc.



Rabbit, dog, ferret plethysmograph



NHP(non-human primate) plethysmograph

## Model Selection


**WBP-4/8/16-MR**

M/R: Mouse/Rat

8: Number of channels, 4/8/16 etc.

# Animal Whole-body plethysmography System


Inhalation toxicology  
studies



Combined with optogenetic  
technology, EEG, EMG, etc.



Simultaneous monitoring of  
animal activity and respiratory



Infrared thermal imaging  
for shell temperature

Comprehensive




Solutions



Hypoxic/CO<sub>2</sub>  
Ventilatory Response



Synchronized video recording



Telemetry system:  
Core temperature, ECG, EEG, etc.



Continuous drug administration,  
blood collection and microdialysis

\*Customization is accepted according to customer experiment requirements

## Publications

- Sun Lingna, Fan Mingrui, Huang Dong, Li Bingqin, Xu Ruoting, Gao Feng, Chen Yanzuo. "Clodronate-loaded liposomal and fibroblast-derived exosomal hybrid system for enhanced drug delivery to pulmonary fibrosis" [J]. *Biomaterials*, 2021
- Jiang C, Huang H, Yang X, et al. Targeting mitochondrial dynamics of morphine-responsive dopaminergic neurons ameliorates opiate withdrawal[J]. *The Journal of Clinical Investigation*, 2024.
  - Jinru Zhang, Keqin Ji, Yuanmeng Ning, Lingna Sun, Mingrui Fan, Chunjie Shu, Ziqi Zhang, Tianyu Tu, Jingyun Cao, Feng Gao, and Yanzuo Chen. Biological Hyperthermia-Inducing Nanoparticles for Specific Remodeling of the Extracellular Matrix Microenvironment Enhance Pro-Apoptotic Therapy in Fibrosis. [J] *ACS Nano* 2023 17 (11)
- Li H, Liu S, Dai W, et al. Pressure-sensitive multivesicular liposomes as a smart drug-delivery system for high-altitude pulmonary edema[J]. *Journal of Controlled Release*, 2024, 365: 301-316.
- Xiong Jiaying, Zhuang Tao, Ma Yurong, Xu Junyi, Ye Jiaqi, Ma Ru, Zhang Shuang, Liu Xin, Liu Bi-Feng, Hao Chao, Zhang Guisen, Chen Yin. "Optimization of bifunctional piperidinamide derivatives as  $\sigma$ 1R Antagonists/MOR agonists for treating neuropathic pain" [J]. *European Journal of Medicinal Chemistry*, 2021, 226:
- Zhang M, Xu B, Li N, et al. All-Hydrocarbon Stapled Peptide Multifunctional Agonists at Opioid and Neuropeptide FF Receptors: Highly Potent, Long-Lasting Brain Permeant Analgesics with Diminished Side Effects[J]. *Journal of Medicinal Chemistry*, 2023.
- Zhang Y, Jiang M, Xiong Y, et al. Integrated analysis of ATAC-seq and RNA-seq unveils the role of ferroptosis in PM2. 5-induced asthma exacerbation[J]. *International Immunopharmacology*, 2023, 125: 111209.
  - Yumei Zhou, Tieshan Wang, Xiaoshan Zhao, Ji Wang and Qi Wang, "Plasma Metabolites and Gut Microbiota Are Associated With T cell Imbalance in BALB/c Model of Eosinophilic Asthma" [J]. *Frontiers*, 2022. 819747:
- Zheng R, Yang W, Wen Yet al. "Dnah9 mutant mice and organoid models recapitulate the clinical features of patients with PCD and provide an excellent platform for drug screening" [J]. *Cell Death Dis* 13, 559 (2022)
- Zhang X, Hu T, Yu X, et al. Human umbilical cord mesenchymal stem cells improve lung function in chronic obstructive pulmonary disease rat model through regulating lung microbiota[J]. *Stem Cells*, 2024: sxae007.
- Wang W, Mu M, Zou Y, Li B, Cao H, Hu D, Tao X. "Inflammation and fibrosis in the coal dust-exposed lung described by confocal Raman spectroscopy" [J]. *PeerJ* 10:e13632
- Dongyi Zhu, Qian Zhang, Qinchuan Li, Guangxue Wang, Zhongliang Guo. "Inhibition of AHNK nucleoprotein 2 alleviates pulmonary fibrosis by downregulating the TGF- $\beta$ 1/Smad3 signaling pathway" [J]. *The Journal of Gene Medicine*, Volume24/Issue9/September 2022/e3442

\*Contact us for more information.