

KN-297 Laboratory animal's ID marker

Opppen



Can be used on the tail and skin of mice and rats

- Available in five colors.
- Capable of marking the tail and skin of not only albino mice but also colored mice such as C57BL/6.
- Animal-safe inks are used.

The method of marking the tail with Opppen is suitable for short-term identification. Compared to ear-punching or ear-marking, which involves cutting the animal's ears, this method is less stressful on the animal and easier for the experimenter.

Twin marker



THICK

THIN


Cat No.	Product Name	for Albino Animals	for Colored Animals
KN-297-B-6	Opppen Black [6pcs/pac]	○	○
KN-297-G-6	Opppen Green [6pcs/pac]	○	○
KN-297-P-6	Opppen Pink [6pcs/pac]	○	○
KN-297-LB-6	Opppen Light Blue [6pcs/pac]	○	○
KN-297-V-6	Opppen Violet [6pcs/pac]	○	○

INTELLIBIO S.A..R.L.
17 avenue du Général de
Gaulle SEICHAMPS -
FRANCE

Tél. 00 + 33 3 83 28 16 04
nelly.parisot@intelli-bio.com
contact@intelli-bio.com

Visitez notre site internet
www.intelli-bio.com

ICR



※The frequency of re-marking varies depending on the strain, individual differences, and age of the mice, etc.
※Color appearance may vary depending on the type of lighting.

Can also be used for marking skin and ears when performing surgery.

Suitable for identification by staining their hair.

INTELLIBIO S.A..R.L.
17 avenue du Général de Gaulle SEICHAMPS - FRANCE

Tél. 00 + 33 3 83 28 16 04
nelly.parisot@intelli-bio.com
contact@intelli-bio.com

Visitez notre site internet
www.intelli-bio.com

